

TOPI-CLICK PERL

Vaginal Application

1) Opening/Attaching the Perl® Snap-On Applicator™

- Flip open the Perl shell cap. (Figure a.)
- Snap on the applicator to the **receptor port** of the Perl dispenser. The applicator plunger should touch the **receptor port**. **The Topi-CLICK Perl dispenser is designed to be used exclusively with the Perl Snap-on Applicator™.** (Figure b.)



2) Loading Medication into the Perl® Snap-On Applicator™

- Simply turn/CLICK the base of the Perl to load the metered dose into the applicator. (Figure c.)
- Keep the Perl in vertical position when loading the applicator as shown. (Figure c.)
- Do not push down on plunger.



Note: Each turn/CLICK dispenses 0.25 mL of medication accurately and does not require fill lines.

Dosing Guide:

1 CLICK=0.25 mL/3 CLICKS=0.75mL/2 CLICKS=0.50mL/4 CLICKS= 1.00mL



3) Removing Perl® Snap-On Applicator™

- Slowly tilt the Perl applicator as shown until the applicator releases from the receptor port. **Do not pull or snap off.** Use vaginal applicator to dose as prescribed. **Applicator is dishwasher safe, but use caution to avoid residual medicine transfer when cleaning.**

4) Closing the Cap

- Flip down shell cap tightly after each dose to prevent evaporation.



Additional Instructions:

- Vaginal cream is best to be used when you are least likely to be intimate within the next 10 hours.
- Avoid sexual intercourse for 10-12 hours after application to prevent transfer of hormone to your partner.
- Do not worry about leakage or 'wasting' the hormone through urination or wiping. The vaginal wall acts like a 'sponge' & will absorb the hormone within the cream.
- Always follow your prescriber's instructions for application when indicated.
- Best used at night so can be absorbed while laying down.
- If using in the morning, consider lying down for 20-30 minutes after so cream is absorbed before standing.