

# Unodose™ Topical Application for Women



1. **Remove** the protective cover with a light counterclockwise twist and pull while holding the barrel.



2. **Turn** the dosing knob in the direction marked for the number of prescribed clicks.



3. **Confirm** that cream has exited the holes of the applicator cap.



4. **Apply** the cream to the skin absorption site as listed on your prescription label. **May Rub** applicator or apply as directed per physician or pharmacist. Rotate areas with each application.



5. **Replace** the protective cover with a light push and clockwise turn. You will hear and feel a click.

*Be sure to leave the knob in a “lock” position, and not between click positions.*



**Check** the scale at the blue line for the number of clicks remaining.

**Refill** your prescription as instructed by your pharmacy.

- Store in a dry place at room temperature.
- Avoid exposure to sunlight and other bright lights.
- When traveling in an airplane or to a different elevation, seal the applicator in a zip lock bag to prevent soiling from any leakage due to pressure changes.