


Unodose™ Vaginal Application



1. **Remove** the protective cover with a light counterclockwise twist and pull while holding the barrel.



2. **Turn** the dosing knob in the direction marked for the number of prescribed clicks.



3. **Confirm** that cream has exited the holes of the applicator cap.



4. **Swipe** cream from top of device onto fingertip. Gently insert fingertip just inside vaginal canal, approximately ½” – 1” depth. Rubbing gently, apply cream in a swiping motion to vaginal.



5. **Replace** the protective cover. *Be sure to leave the knob in a “lock” position, and not between click positions.* **Check** the scale at the blue line for the number of clicks remaining. **Refill** your prescription as instructed by your pharmacy. Store in a dry place at room temperature.

Precautions

- 1) **Avoid sexual intercourse for 10-12 hours after application to prevent transference of hormone to your partner.**
- 2) **The vaginal wall acts like a ‘sponge’ & will absorb the hormone within the cream. Therefore, you need not worry about ‘wasting’ the hormone through urination, wiping with tissue paper, or leakage.**
- 3) **Always follow your prescriber’s instructions for application when indicated.**

Avoid exposure to sunlight and other bright lights. When traveling in an airplane or to a different elevation, seal the applicator in a zip lock bag to prevent soiling from any leakage due to pressure changes.